

EXAMINATION LAMP

LP01-28W



Lampu *LED*



Daya Lampu 28 W



Lifetime 50.000 jam



Desain yang ekonomis



Pengoperasian yang mudah

TKDN
48.49%



Scan Me!

Nama Produk	Examination Lamp
Tipe	LP01-28W
Kelas	1
IP	IPX1
Tegangan input	110 - 240 VAC 50/60 Hz
Daya lampu <i>LED</i>	28 W
Temperatur warna	6000±500K
<i>CRI</i>	80%
Suhu lingkungan	5°C - 40°C
Kelembapan	≤85%
Tekanan atmosfer	86.0 - 106 kPa
<i>Lifetime</i>	50.000 jam