

Examination Lamp

LP02-9W



Lampu LED



Daya Lampu 9 W



Desain yang ekonomis



Pengoperasian yang mudah



Nama Produk	MIKI Examination Lamp
Tipe Produk	LP02-9W
Kelas Perlindungan Produk	1
Tingkat IP	IPX1
Daya AC	110~220 VAC, 50/60 Hz
Daya Lampu	9 W
Temperatur Warna	6000 $\pm$ 500 K
Suhu Lingkungan	5 °C~40 °C
Kelembaban Relatif	$\leq$ 85%
Tekanan Atmosfer	86.0~106 kPa