

SUCTION PUMP SP-05S

TKDN
42.45%



Dapat menghisap cairan tubuh seperti darah, air liur, lendir, dan lain-lain



Mesin Vakum tanpa minyak dan pelindung termal, *auto reset*



Tabung Suction berbentuk kapsul yang dapat mencegah luapan cairan



Catheter



Filter



Silicone Tube

KEMENKES RI AKD : 21603422861



ISO 9001:2015



ISO 13485:2016
ISO 45001:2018
ISO 14001:2015



Scan Me!



SPESIFIKASI

Nama alat	MIKI Suction Pump
Tipe alat	SP-05S
Rentang vakum	>650 mmHg
Laju aliran maksimum	40 liter/menit
Kapasitas	1500 c.c.
Pengaman alat	Kontrol Pengambang
Meter kebisingan	60 - 65 dB
Daya	220 V
Aliran	1.1 A / 0.85 A
R.P.M	1720 rpm
Keluaran	80 W
Ukuran alat	34 x 17 x 40 cm
Ukuran <i>packing</i>	37 x 24 x 35.5 cm
Berat	5 kg

Diproduksi Oleh :
PT MENARA INTI KREASI INDONESIA
Tangerang - Banten
 mikiindo.com

Didistribusi Oleh :